

**Table 3:** Carries contamination rate results regarding the chemiluminescence analyzer performance based on the Clinical and Laboratory Standards Institute H57-A document requirements.

<b>Test items</b>	<b>L1-L3</b>	<b>H3-L3</b>	<b>Carried contamination rate (%)</b>
TM	0.20	165.40	0.12
TAT	-0.10	92.30	-0.11
PIC	0.02	19.40	0.10
t-PAIC	0.20	46.70	0.43

TM: thrombomodulin, TAT: thrombin-antithrombin complex, PIC: plasmin- $\alpha$ 2-plasmin inhibitor complex, tPAIC: tissue plasminogen activator/plasminogen activator inhibitor-1.